The Super B/W Laptop Ultrasound System

CTS-900



CTS-900 is

The upgrade of SIUI laptop B/W ultrasound product



New generational chip and advanced signal processing technology

- Strongest integrate circuit, highest capability of transmit and receipt
 - →good performance and steady stability
- New apparatus and components
 - → compact hardware and exquisite system (main unit with battery:2.8kg!!!)



New generational chip and advanced signal processing technology

- original creation: segment logarithmic compression →smooth image with less grain
- AD convert with high velocity and resolution
 - → excellent penetration
- Dynamic frequency transmit technology
 - → deep scan depth up to 25.2cm, suitable for clinical diagnosis



Characteristics and Features

- 10.4 inch LCD monitor with high resolution, excellent contrast, anti-dazzling
 - →ensure good vision in any environment.
- Built-in 2G CF card, 2 USB ports and temporary store
 - →strong storage capacity
- One key optimization
 - → User-defined, 14 kinds of parameter adjustable
- Optional DICOM3.0
 - → Image Transmit



Characteristics and Features

Power saving function

Probes auto freeze

→ Extend service life of probe effectively

Screen saver

→ Extend service life of the LCD monitor

System hibernation

→ Save battery power



Characteristics and Features

- The system supports more than 11 kinds of broadband probes including convex, linear, micro-convex, vaginal, rectal probes
 - →meet great diagnosis demand of different organs and tissues
 - → clear and high resolution images for superficial diagnosis



Characteristics and Features

Complete application software

- A. Obstetrics
- B. Gynecology
- C. Cardiology
- D. Small parts
- E. Abdomen

- F. Urology
- G. Pediatrics
- H. Orthopedics
- I. Peripheral Vessel
- J. Podiatry




the premium 2D imaging performance - liver





the premium 2D imaging performance - Liver





the premium 2D imaging performance – Kidney






the premium 2D imaging performance -Thyroid






the premium 2D Chroma imaging performance

→ 9 kinds for selection






the premium 2B imaging performance





the premium 4B imaging performance





Special Standard configuration

- Battery
 - → at less 1.5 hours up to 2 hours (depend on different application modes)
- Power Supply Adapter
 - →better stability
- Car charger
 - →convenient to charge on driving










Particular Optional Accessories



CTS-900 –SIUI B/W Laptop Ultrasound Particular Optional Accessories

Carry case → protection






Layout of the carry case



→ Main unit

Lower Module

- → Adapter
- → Car charge
- → Cables
- → Accessory






Layout of the lower module



battery

CTS-900 -SIUI B/W Laptop Ultrasound

Layout of the upper module





Contact us for further information



www.siui.com



siui@siui.com



86-754-88251499



86-754-88250150

